

# Mobility Solutions



Blue  
Gray  
Rose

#2000

Since 1986, the Columbia Orthopedic Positioning Seat has been parents' and therapists' first choice in special needs car seats. Appropriate for children 20 to 102 lbs and up to 5' tall.

- **Lightest car seat of its size** – only 14 lbs.  
Easy to use with optional mobility bases (see backside).
- **Toughest** Meets or exceeds safety standards: American FMVSS 213, Canadian CMVS 213.3 and European ECE R44.03. FAA compliant for air travel.
- **Convenient** Easy to install; fits all cars, vans, school buses and most airplane coach seats.
- **Most Comfortable** Luxurious –and durable–padded, machine-washable cover, in 3 attractive colors.
- **Grows with the child!**
  - Height adjustable head supports and 5-point harness
  - Seat extensions add depth of 2 1/2" and 5"
  - Extensor Thrust Wedge changes seat to back angle
  - Lateral positioning pads narrow the seat width from 15" down to 8"
  - Block style pads for special positioning needs



## Included Accessories

- 5-Point Harness • Shoulder Harness Pads
- Lateral Positioning Pads ...custom fit and comfort for the small or slender child
- Adjustable Head Pads • Crotch Strap Pad • Tether Strap Kit
- Easy-to-follow Installation Video & Instructions, in English with Spanish subtitles

## Additional Accessories



Closed

#2050, #2055 Wedge



changes seat-to-back angle



#2025 Seat Depth Extender adds 2 1/2"



#2026 Seat Depth Extender adds 5"



Open

#2015 Swing-Away Abductors in two heights for custom comfort



#2075 Buckle Guard keeps "little Houdini" buckled

#2008 4 Block Style Positioning Pads  
Two pads 2.5" thick, Two pads 1.5" thick  
(Specify color)



Mobility Solutions  
from Columbia Medical

#2175 Handy Tote attaches to back rail of Base



### Tilt-in-Space Mobility Base # 2100

The Tilt-in-Space Mobility Base is a strong, lightweight frame custom designed for the Columbia Orthopedic Positioning Seat. Perfect for shopping at the mall, going to the playground or when traveling.

- **Comfortable**—five angles of tilt and height-adjustable footrest.
- **Easy to use**—seat locks easily onto the base.
- **Convenient**—folds to carry, or store compactly.
- **Durable**—strong yet lightweight frame.
- **Built to go everywhere**—large wheels for a smooth ride over any surface.

### Special Needs Baby Jogger # 2200



The Special Needs Baby Jogger can be used alone or with the Columbia Orthopedic Positioning Seat. Designed with convenience and maneuverability in mind, the Special Needs Baby Jogger opens up a new world for the physically challenged. Its bike-style wheels roll easily over gravel, sand, trails and other terrain normally inaccessible to conventional wheelchairs. For city strolling, its front wheel swaps for swivels to maneuver through crowds and shops. Standard features include a retractable sun canopy, handbrake, footbrake, five-point safety harness, and footplate.


### RunAbout by Convoid # 2300

Go anywhere with the RunAbout by Convoid. Use alone or with the Columbia Orthopedic Positioning Seat. Compact, solid and durable, the lightweight portability makes it a convenient answer for busy mothers and children.

- Accessories:** Heel Loops (#2301)  
Ankle Cuffs (#2302)  
5-point H-Harness (#2303)  
Canopy (#2304)



## Product Specifications



### Orthopedic Positioning Car Seat in Blue, Gray, Rose

<b>Model 2000</b>	
User weight	20-102 lbs (9-46kg)
User height	up to 5' (152cm) tall
Product weight	14 lbs (6kg)
Back height	29" (74cm)
Seat (inside)	15"W x 10"D (38x25cm)


### Positioning Pads

4 Lateral pads included with seat reduce seat width from 15" to 12"

<b>Model 2010 -Lateral style</b>	
Extra set of 4 pads	
<b>Model 2008 -Block style</b> specify upholstery color	
Set of 4 pads (two 2.5" thick and two 1.5" thick)	

### Seat Wedge

<b>Model 2050</b>	
2" (5cm) front height slopes to 1" (3cm) rear	
<b>Model 2055</b>	
3" (8cm) front height slopes to 1" (3cm) rear	



### Growth Adjustments

Abductor	Seat Depth	Seat Extender
Model: 2015 or 2015T	10" (25cm)	None
Model: 2018 or 2018T	12.5" (31cm)	Model: 2025
Model: 2019 or 2019T	15" (39cm)	Model: 2026

### Seat Depth Extender

<b>Model 2025</b> increases seat depth to 12.5" (30cm)
<b>Model 2026</b> increases seat depth to 15" (38cm)

### Swing-Away Abductor: Standard & Tall

<b>Model 2015, 2018, 2019</b>	3.5"W x 4"H x 3"D
<b>Model 2015T, 2018T, 2019T</b>	3.5"W x 6"H x 3"D

### Tilt-in-Space Mobility Base

<b>Model 2100</b>	
Five angles of tilt	20° 25° 30° 35° 40° from vertical
Open	22.5"W x 47"L x 37"H(57x119x94cm)
Folded	22.5"W x 43"L x 19"H(57x109x48cm)
Leg drop	5"-16"(13-41cm)
Front wheels	7"(18cm) diameter
Rear wheels	11"(28cm) diameter
Product weight	25 lbs(11kg)
Seat width	16"(41cm)



800-454-6612 ext. 100

email: [info@columbiamedical.com](mailto:info@columbiamedical.com) • [www.columbiamedical.com](http://www.columbiamedical.com)